

## JC

### YELLOW CROWN AND BRIDGE ALLOY

JC is a type III crown and bridge casting gold. It is one of a family of seven high gold content casting alloys offered by Jensen Industries. JC features 78% nobility, rich gold color, high strength, competitive cost, and melts, casts, finishes and polishes with speed, ease and accuracy.

#### PROPERTIES

Melting Range .....	1670 to 1760°F
Density .....	15.2 g/cm <sup>3</sup>
Grain Size .....	24 microns
Hardness .....	125 HV
Tensile Elongation .....	50%
Tensile Yield Strength (PSI) .....	39,500
Ultimate Tensile Strength (PSI) .....	62,000

#### CHEMISTRY

Gold .....	75%
Silver .....	11%
Copper .....	10%
Palladium .....	3 %

Contains less than 1%  
Zinc, Indium, Iridium

Classification - High Noble

#### PROCESSING TECHNIQUE

##### SPRUNG

The indirect method is recommended for multi-units. Use an 8 gauge runner bar with 10 gauge connectors. If preferred, the direct method may be used on both single units and small bridges. Use a 10 gauge sprue 1/4" (6mm) to 3/8" (9mm) long. Sprues longer than 3/8" (9mm) should have a reservoir 1/16" (1.5mm) from pattern. Patterns should be a maximum of 1/4" (6mm) from top of investment.

##### INVESTMENT AND BURNOUT

Either gypsum or phosphate bonded investment may be used following the manufacturer's instructions. The burnout temperature should be at least 900°F (480°C) and should not exceed 1200°F (650°C).

##### MELTING AND CASTING

Extra winds of the casting arm are not required. A gas/compressed air or gas/oxygen flame with 5 PSI gas and 10 PSI oxygen is recommended. The alloy will fully puddle and form a ball before it is ready to cast. The casting temperature is 1850°F (1010°C). Bench cool to a dull red heat before quenching.

##### DEVESTING AND FINISHING

Blast with aluminum oxide to remove investment particles and oxidation. Finish and polish using standard techniques.

##### SOLDER AND FLUX

Solder:	650 Fine Solder
Flux:	Brown Fluoride Flux